



The future of drug screening is automated with the RapidReader App


Streamline your drug screening process by alleviating the burden of physical paperwork and manual inputs. The RapidReader App is a reliable solution to support users with results collection, reporting, and upload directly from any mobile device.



Streamline your drug screening workflow with the power of AI.

The **RapidReader App** is a state-of-the-art screening software application that allows users to streamline their drug screening process. The RapidReader App uses AI technology to capture drug screen results from a picture, record the results and generate a Chain of Custody form.

The **RapidReader App** is compatible with a wide range of our Rapid Response™ drug screening products - trusted for over a decade to provide fast, reliable and accurate results.



Efficient

- Accelerated workflow
- Average processing time of just 2½ minutes
- Accessible from anywhere



Automated

- Innovative AI technology
- Capture & record results directly from test cups and panels



Streamlined

- Intuitive user experience
- Securely review & confirm results
- Integrates seamlessly with users' Internal Data Management Systems (ATS/EMR)

Collect a touchless signature on the spot

Save time by having participants sign from their own devices using our touchless signature feature.

< Back Consent Form

Donor Info
ID: 560112
Name: Rapid Response
Date-of-Birth: 7 Oct 2001

I authorize BTNX Inc. to perform a test.

Use Touchless Signature: ☐ [clear](#)

Rapid Response
sign above the line

Request Signature Continue with test

Capture and record results from a photo

Take a photo of the results window on the drug screening cup to quickly capture results for review.

< Drug Screening - Saliva Test

Roger Response

Specimen ID: SZ021070043

Rapid Response
DIAGNOSTICS

Saliva Screening Test
ID: _____
Date: _____

CAPTURE

Review and confirm results in seconds

Using the data from the photo, the built-in AI software autofills results for review, confirmation, and update.

Drug Screen Test Results

Your Company

COC	Negative	Edit
MET	Negative	Edit
AMP	Negative	Edit
OPI	Non-Negative	Edit
OXY	Negative	Edit
MTD	Negative	Edit
PCP	Negative	Edit
BAR	Non-Negative	Edit
BZO	Negative	Edit
MDMA	Negative	Edit

Notes:

Notes

IMAGES

ALCOHOL TEST

CONTINUE

Share and upload paperless Chain of Custody forms

Results, signatures, donor information and test images can be emailed to stakeholders or uploaded to your organization's internal data management system, ATS, or EMR.

BTNX inc.

Rapid Response

DRUG SCREEN

DRUG AND ALCOHOL SCREEN

Chain of Custody Form

CLIENT NAME: Jane Do
570 Hood Rd. Markham ON L3R 4G7

OPERATOR: Roger Response

CANDIDATE: 1007

DRUG SCREEN

☒ Urine Test ☐ Saliva Test

Lot #: DOA2109599 Lot Expiry: 9/24/2023 Test Date: 2022-10-14 Reason for Screen: Pre employment

DRUG	RESULTS	DRUG	RESULTS
COC	NEGATIVE	PCP	NEGATIVE
MET	NEGATIVE	BAR	NON-NEGATIVE
AMP	NEGATIVE	BZO	NEGATIVE
OPI	NON-NEGATIVE	MDMA	NEGATIVE
OXY	NEGATIVE		
MTD	NEGATIVE		

DRUG TEST RESULTS:

Positive

This screen has produced a Positive result for one or more of the parameters listed above.

ALCOHOL SCREEN

Lot #: Lot Expiry: Test Date: Reason for Screen:

ALCOHOL TEST RESULTS:

☐ Negative Screen ☐ Positive Screen

Level: ☐ 0% ☐ 0.2% ☐ 0.4% ☐ 0.8% ☐ 3%

EMPLOYEE SIGNATURE: Roger Response

OPERATOR SIGNATURE: Jane Do

WITNESS SIGNATURE:

This document contains personal identifiable information that must be treated confidentially. Any unauthorized use or disclosure is prohibited.

Page 1 of 2

www.btnx.com

BTNX inc.

Rapid Response

DRUG SCREEN

DRUG AND ALCOHOL SCREEN

Chain of Custody Form

CLIENT NAME: Jane Do
570 Hood Rd. Markham ON L3R 4G7

OPERATOR: Roger Response

CANDIDATE: 1007

DRUG SCREEN TEST IMAGE(S)

Lot #: DOA2106359 Lot Expiry: 9/24/2023 Test Date: 2022-10-14 Reason for Screen: Pre employment

BTNX inc.

Rapid Response

DRUG SCREEN



“ The RapidReader App has been a game changer for our business. In a world of evolution and rapidly advancing technology along with the high volume we experience, it was imperative we move to an automated drug screen process and that’s exactly what we did with Lochness. With the Rapid Reader App, **we were able to drastically decrease the amount of time it takes to onboard while increasing our candidate and recruiter experience.** They also provide excellent customer service and are always available when we need them. I highly recommend Lochness for all of your screening needs!

– Hailey Fernandez, Director of Corporate Services with Effex Management Solutions ”

5X Faster!

Process test results and generate a Chain of Custody in just 2 ½ minutes, 5x faster than manual inputs*. Save time and streamline your screening process with automated results.

*Based on time to record results and create equivalent paper chain of custody documents.



Using an ATS or EMR system?

We offer custom integrations with ATS and EMR systems, allowing your team to save screening results to your system with the click of a button, no complicated file transfers required.



Onboarding is a breeze

We offer comprehensive training resources, demos, and ongoing technical support to ensure every member of your team can screen with confidence.

Rapid Response™ Saliva Cup

The Rapid Response™ Saliva Test Cup is a rapid, non-invasive drug screening solution that can simultaneously detect multiple drugs in saliva samples.

Specimen: **Saliva Sample**

- Results in just 10 minutes
- Available in a wide range of drug parameter combinations
- Easy to monitor collection to prevent adulteration
- The RapidReader App is 100% compatible with Rapid Response™ brand Saliva Test Cups



Rapid Response™ Multi-Drug Urine Tests

The Rapid Response™ Flat Test Cup and Rapid Response™ Test Panel are easy-to-use drug screening solutions that can simultaneously detect multiple drugs in urine samples.

Specimen: **Urine**

- Results in just 5 minutes
- Available in a wide range of drug parameter combinations
- The RapidReader App is 100% compatible with Rapid Response™ brand Flat Test Cups and Test Panels

Flat Test Cup

- No direct contact with sample
- Leak-proof cap
- Built-in temperature strip

Panels

- Compact format
- Easy to store and transport





**Ready to streamline your
drug screening process with
the RapidReader App?**



**Scan to
see the app
in action**

Equipping your staff with an
all-in-one screening solution has never been easier.
Our team at **Lochness Medical Supplies** is ready to
walk you through the process. Book a demo with us at:

Tel: +1 888-506-2658

Email: info@lochnessmedical.com

www.lochnessmedical.com



© 2023 BTNX Inc.

All rights reserved. Rapid Response™ is a registered
trademark of BTNX Inc.